

**Improved liquid chromatography-tandem mass spectrometric method for the determination of ethyl glucuronide concentrations in hair: Applications to forensic cases**

Imbert L, Gaulier JM, Dulaurent S, Morichon J, Bevalot F, Izac P, Lachâtre G.

International journal of legal medicine

2014; 128(1):53-58

**ARTICLE IDENTIFIERS**

DOI: 10.1007/s00414-013-0894-8

PMID: 23824336

PMCID: not available

**JOURNAL IDENTIFIERS**

LCCN: not available

pISSN: 0937-9827

eISSN: 1437-1596

OCLC ID: 23022400

CONS ID: sn 91033072

US National Library of Medicine ID: 9101456

This article was identified from a query of the SafetyLit database.